



CORFLEX

**CORFLEX IS A GLOBALLY BASED
MANUFACTURER AND DISTRIBUTOR OF QUALITY
ORTHOPEDIC AND REHABILITATION PRODUCTS**



SHOULDER IMMOBILIZERS

ULTRA SHOULDER IMMOBILIZER

- > Features 1½" wide waist strap for secure immobilization and thumb loop to help prevent wrist drop
- > Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery, or other shoulder related injuries
- > Manufactured from seamless, lightweight, durable material with added shoulder strap pad



PRODUCT NUMBER

2XS	XS	S	M	L	XL
23-1869	23-1870	23-1871	23-1872	23-1873	23-1874

SIZING GUIDE

	2XS	XS	S	M	L	XL
Width/Height	10" x 4"	12" x 5"	15" x 6"	17" x 7"	19" x 8"	21" x 9"

Measurements are sling pouch width x height

APPLICATION INSTRUCTIONS

1. Open sling and place arm in inside, making sure elbow is in sling as far back as possible.
2. Thread strap through d-rings.
3. Secure hook and loop so elbow is positioned at 90-degree angle. (If necessary, there is a one-time hook and loop adjustment behind the elbow to lengthen or shorten strap length as needed.)
4. Position elbow (shoulder) strap and pull behind your back.
5. Thread strap through buckle and pull tight to prevent arm from moving away from body.

CARE & CLEANING

Corflex manufactures a quality line of products and to help your care for the product properly we provided the following instructions for care & cleaning:

Close all contact closures and hand wash in warm water with mild soap. Product should be air-dried, do not use high temperatures to wash or dry.

HOW TO PLACE AN ORDER

The Corflex Customer Service Team is ready to assist you in placing Your order. You may place your orders via a variety of methods:



TELEPHONE

1-800-426-7353, Ext1



FAX

1-603-623-4111



E-MAIL

Service@corflex.com



MAIL

CORFLEX, Inc.
Attn: Customer Service
669 East Industrial Park Dr
Manchester, NH 03109 USA

669 East Industrial Park Drive - Manchester, NH 03109
Phone (603) 623-3344 - Fax (603) 623-4111 - www.corflex.com

© Corflex | The advice and information contained in this document is not intended as a substitute for medical counseling. If you have any questions about treatment for a specific condition, please check with a health care professional.