



CORFLEX

**CORFLEX IS A GLOBALLY BASED
MANUFACTURER AND DISTRIBUTOR OF QUALITY
ORTHOPEDIC AND REHABILITATION PRODUCTS**



SHOULDER IMMOBILIZERS

ULTRA ER SHOULDER ABDUCTION PILLOW W/SLING

- > Features removable wedge that provides 30° of external rotation and V-lock strapping system to prevent rolling and shifting
- > Indicated for protocols involving non-surgical immobilization in external rotation to help reduce the incidence of recurrent dislocation
- > Constructed of firm density foam covered in tricot material with tricot sling and exercise ball to stimulate circulation



PRODUCT NUMBER

S	M	L	XL
23-1941	23-1942	23-1943	23-1944

SIZING GUIDE

	S	M	L	XL
Hip	24"-30"	30"-36"	36"-42"	42"-48"
Forearm	9"-12"	12"-15"	15"-18"	18"-20"

Measurement taken from hip circumference and length of forearm from wrist to elbow

APPLICATION INSTRUCTIONS

1. Grasp pillow assembly, making sure wider end of wedge is facing forward, release buckle and place curved side of pillow assembly on waistline under affected shoulder.
2. Grasp both ends of waist strap, secure buckle in front and adjust to desired tightness around waist (Note: to trim excess strap material open alligator tip, trim and replace alligator tip).
3. Remove sling from plastic pouch, make sure opening is facing forward/upward, line up with hook and affix to side of pillow assembly.
4. Open sling, place arm in sling and position elbow in the sling as far back as possible, and place thumb through loop inside sling for additional comfort.
5. Grasp sling strap, bring around back and place unaffected arm through opening so strap is comfortably placed at base of neck, around shoulder and under arm.
6. Feed end of sling strap through front sling buckles, secure strap and close top of sling to secure around arm.
7. An optional exercise ball w/removable holder is attached to the shoulder abduction pillow wedge and may be used to encourage circulation.

NOTE: The pillow is designed to keep the arm in approx. 30 deg. of abduction. The position of the pillow at the waist should be positioned for desired internal or external rotation by sliding the pillow forward or back along the waistline. Also, wedge portion of pillow may be removed if external rotation is not desired.

CARE & CLEANING

Corflex manufactures a quality line of products and to help your care for the product properly we provided the following instructions for care & cleaning:

Remove any metal stays if possible, close all contact closures and hand wash in warm water with mild soap. Product should be air-dried, do not use high temperatures to wash or dry.

HOW TO PLACE AN ORDER

The Corflex Customer Service Team is ready to assist you in placing Your order. You may place your orders via a variety of methods:



TELEPHONE

1-800-426-7353, Ext1



FAX

1-603-623-4111



E-MAIL

Service@corflex.com



MAIL

CORFLEX, Inc.
Attn: Customer Service
669 East Industrial Park Dr
Manchester, NH 03109 USA

669 East Industrial Park Drive - Manchester, NH 03109
Phone (603) 623-3344 - Fax (603) 623-4111 - www.corflex.com

© Corflex | The advice and information contained in this document is not intended as a substitute for medical counseling. If you have any questions about treatment for a specific condition, please check with a health care professional.