



**CORFLEX**

**CORFLEX IS A GLOBALLY BASED  
MANUFACTURER AND DISTRIBUTOR OF QUALITY  
ORTHOPEDIC AND REHABILITATION PRODUCTS**



## BACK SUPPORTS

### CRISS-CROSS BACK SUPPORT

- > Features criss-cross back panel with flexible anti-roll stays and double-pull elastic tension straps for added compression
- > Indicated for moderate low back pain, lumbar muscle weakness and lumbar sprains/strains
- > Manufactured from premium quality elastic with contact closure



#### PRODUCT NUMBER

XS	S	M	L	XL	2XL
35-1600	35-1601	35-1602	35-1603	35-1604	35-1605
3XL	4XL				
35-1606	35-1607				

#### SIZING GUIDE

XS	S	M	L	XL	2XL
18"-24"	24"-30"	30"-36"	36"-42"	42"-48"	48"-52"
3XL	4XL				
52"-56"	56"-60"				

Measurement taken from hip circumference

### APPLICATION INSTRUCTIONS

1. Lay back support flat, locate top of support by identifying openings for stays (top of support). If applicable, stay assemblies may be removed and adjusted.
2. Grasp back support and center on small of back.
3. Pull left side firmly across middle section, repeat for right side. If applicable, attach double pull and secure.

### CARE & CLEANING

Corflex manufactures a quality line of products and to help your care for the product properly we provided the following instructions for care & cleaning:

Remove any metal stays if possible, close all contact closures and hand wash in warm water with mild soap. Product should be air-dried, do not use high temperatures to wash or dry.

### HOW TO PLACE AN ORDER

The Corflex Customer Service Team is ready to assist you in placing Your order. You may place your orders via a variety of methods:



**TELEPHONE**

1-800-426-7353, Ext1



**FAX**

1-603-623-4111



**E-MAIL**

Service@corflex.com



**MAIL**

CORFLEX, Inc.  
Attn: Customer Service  
669 East Industrial Park Dr  
Manchester, NH 03109 USA

669 East Industrial Park Drive - Manchester, NH 03109  
Phone (603) 623-3344 - Fax (603) 623-4111 - www.corflex.com

© Corflex | The advice and information contained in this document is not intended as a substitute for medical counseling. If you have any questions about treatment for a specific condition, please check with a health care professional.