



**CORFLEX**

**CORFLEX IS A GLOBALLY BASED  
MANUFACTURER AND DISTRIBUTOR OF QUALITY  
ORTHOPEDIC AND REHABILITATION PRODUCTS**



## WRIST SPLINTS

### INDUSTRIAL WRIST SPLINT

- > Low profile design provides flexible support, making it ideal for the work place
- > Indicated for Carpal Tunnel Syndrome, tendonitis, strains or sprains, repetitive stress injuries
- > Manufactured from open-weave material, flexible stays and heavy-duty circumferential elastic wrist strap



#### PRODUCT NUMBER

	XS	S	M	L	XL
8" Right	73-5620	73-5621	73-5622	73-5623	73-5624
8" Left	73-5630	73-5631	73-5632	73-5633	73-5634

#### SIZING GUIDE

XS	S	M	L	XL
4"-5"	5"-6"	6"-7"	7"-8"	8"-9"

Measurement taken from wrist circumference

### APPLICATION INSTRUCTIONS

1. Open wrist splint and lay flat (inside up) on any clean surface.
2. Place hand on splint so that the strap with the loop portion through the index finger and the thumb and attaches to hook on top of your hand.
3. Pull thin strap between index finger and thumb and attach to hook closure on top of splint.
4. Wrap 3" elastic strap around wrist, securing to desired comfort level.

### CARE & CLEANING

Corflex manufactures a quality line of products and to help your care for the product properly we provided the following instructions for care & cleaning:

Remove any metal stays if possible, close all contact closures and hand wash in warm water with mild soap. Product should be air-dried, do not use high temperatures to wash or dry.

### HOW TO PLACE AN ORDER

The Corflex Customer Service Team is ready to assist you in placing Your order. You may place your orders via a variety of methods:



**TELEPHONE**

1-800-426-7353, Ext1



**FAX**

1-603-623-4111



**E-MAIL**

Service@corflex.com



**MAIL**

CORFLEX, Inc.  
Attn: Customer Service  
669 East Industrial Park Dr  
Manchester, NH 03109 USA

669 East Industrial Park Drive - Manchester, NH 03109  
Phone (603) 623-3344 - Fax (603) 623-4111 - www.corflex.com

© Corflex | The advice and information contained in this document is not intended as a substitute for medical counseling. If you have any questions about treatment for a specific condition, please check with a health care professional.