



CORFLEX

**CORFLEX IS A GLOBALLY BASED
MANUFACTURER AND DISTRIBUTOR OF QUALITY
ORTHOPEDIC AND REHABILITATION PRODUCTS**



ANKLE IMMOBILIZERS

LOWER LEG FIXED METAL WALKER

- > Provides stability for lower leg; low natural heel height reduces leg length discrepancy
- > Indicated for stable foot/ankle fractures, severe sprains or strains, soft tissue injuries, post-op and post-cast removal
- > Constructed of rigid malleable metal uprights, padded liner and contact closure straps



PRODUCT NUMBER

XS	S	M	L	XL
79-1049	79-1050	79-1051	79-1052	79-1053

SIZING GUIDE

	MENS SIZE	WOMENS SIZE
X-Small	up to 5	up to 5
Small	5 - 6	5 - 7
Medium	6 - 10	7 - 11
Large	10 - 12	11 - 13 1/2
X-Large	12 +	14 +



Measured by shoe size

APPLICATION INSTRUCTIONS

1. Remove the liner from the boot and apply with the ankle at 90 degree angle and the foot flat on the ground.
2. Pull the uprights apart and slide the leg back into the walker frame. Take care to position the ankle at 90 degrees and the uprights centered up the medial and lateral sides of the lower leg.
3. Starting with the bottom strap fasten all the straps from bottom to top.

CARE & CLEANING

Corflex manufactures a quality line of products and to help your care for the product properly we provided the following instructions for care & cleaning:

Wipe plastic shell with damp cloth, you may use disinfectant or mild soap, and wipe dry.

For products with a removable liner, you may hand wash liner in warm water with mild soap, air dry, do not use high temperatures to wash or dry.

HOW TO PLACE AN ORDER

The Corflex Customer Service Team is ready to assist you in placing Your order. You may place your orders via a variety of methods:



TELEPHONE

1-800-426-7353, Ext1



FAX

1-603-623-4111



E-MAIL

Service@corflex.com



MAIL

CORFLEX, Inc.
Attn: Customer Service
669 East Industrial Park Dr
Manchester, NH 03109 USA

669 East Industrial Park Drive - Manchester, NH 03109
Phone (603) 623-3344 - Fax (603) 623-4111 - www.corflex.com

© Corflex | The advice and information contained in this document is not intended as a substitute for medical counseling. If you have any questions about treatment for a specific condition, please check with a health care professional.