



**CORFLEX**


**CORFLEX IS A GLOBALLY BASED  
MANUFACTURER AND DISTRIBUTOR OF QUALITY  
ORTHOPEDIC AND REHABILITATION PRODUCTS**



## FINGER SPLINTS

### STAX SPLINTS

- > Designed to provide immobilization of the distal finger or thumb while allowing for examination of finger tip and nailbed
- > Indicated for treatment of finger tip injuries or dislocations, nailbed injuries, mallet finger and protection following injury
- > Manufactured from durable plastic



#### PRODUCT NUMBER

#1	#2	#3	#4	#5	#5.5
49-3801	49-3802	49-3803	49-3804	49-3805	49-3806

#6	#7	KIT
49-3807	49-3808	49-3800

#### SIZING GUIDE

#1	#2	#3	#4	#5	#5.5
2"	2 1/8"	2 1/4"	2 1/2"	2 5/8"	2 3/4"

#6	#7	KIT
2 7/8"	3 1/4"	30pcs Assortment

Measurement taken at DIP joint circumference

### APPLICATION INSTRUCTIONS

1. Select correct size.
2. Contour aluminum to individual patient based on condition.
3. Secure splint to hand.

### CARE & CLEANING

Corflex manufactures a quality line of products and to help your care for the product properly we provided the following instructions for care & cleaning:

Wipe plastic shell with damp cloth, you may use disinfectant or mild soap, and wipe dry.

For products with a removable liner, you may hand wash liner in warm water with mild soap, air dry, do not use high temperatures to wash or dry.

### HOW TO PLACE AN ORDER

The Corflex Customer Service Team is ready to assist you in placing Your order. You may place your orders via a variety of methods:



**TELEPHONE**

1-800-426-7353, Ext1



**FAX**

1-603-623-4111



**E-MAIL**

Service@corflex.com



**MAIL**

CORFLEX, Inc.  
Attn: Customer Service  
669 East Industrial Park Dr  
Manchester, NH 03109 USA

669 East Industrial Park Drive - Manchester, NH 03109  
Phone (603) 623-3344 - Fax (603) 623-4111 - www.corflex.com

© Corflex | The advice and information contained in this document is not intended as a substitute for medical counseling. If you have any questions about treatment for a specific condition, please check with a health care professional.