



CORFLEX

**CORFLEX IS A GLOBALLY BASED
MANUFACTURER AND DISTRIBUTOR OF QUALITY
ORTHOPEDIC AND REHABILITATION PRODUCTS**



FOOT / ANKLE MISC

RIGID POST-OP SHOE

- > Features flat footbed with rigid sole to prevent excessive foot motion while walking
- > Indicated for post-op foot protection and support
- > Constructed of wood-like sole with water-resistant upper and contact closure



PRODUCT NUMBER

	PEDI	S	M	L	XL
Mens	69-3122	69-3123	69-3124	69-3125	
Womens	69-3132	69-3133	69-3134		
Pedi	69-3121				

SIZING GUIDE

LENGTH	PEDI	S	M	L	XL
Mens		11"	12"	13"	14"
Womens		9"	10"	11"	
Pedi	8"				



APPLICATION INSTRUCTIONS

1. Place affected foot inside shoe.
2. Wrap the straps over the product onto opposite side, and attach straps onto hook portion of the product.
3. Straps should be pulled gently for comfortable fit.

CARE & CLEANING

Corflex manufactures a quality line of products and to help your care for the product properly we provided the following instructions for care & cleaning:

Close all contact closures and hand wash in warm water with mild soap. Product should be air-dried, do not use high temperatures to wash or dry.

HOW TO PLACE AN ORDER

The Corflex Customer Service Team is ready to assist you in placing Your order. You may place your orders via a variety of methods:



TELEPHONE

1-800-426-7353, Ext1



FAX

1-603-623-4111



E-MAIL

Service@corflex.com



MAIL

CORFLEX, Inc.
Attn: Customer Service
669 East Industrial Park Dr
Manchester, NH 03109 USA

669 East Industrial Park Drive - Manchester, NH 03109
Phone (603) 623-3344 - Fax (603) 623-4111 - www.corflex.com

© Corflex | The advice and information contained in this document is not intended as a substitute for medical counseling. If you have any questions about treatment for a specific condition, please check with a health care professional.