



CORFLEX

**CORFLEX IS A GLOBALLY BASED
MANUFACTURER AND DISTRIBUTOR OF QUALITY
ORTHOPEDIC AND REHABILITATION PRODUCTS**



THUMB SPLINTS

THUMBSTER SOFT

- > Anatomically designed to provide comfortable and secure immobilization of the thumb and wrist support while preserving mobility of the other fingers
- > Indicated for deQuervain's Syndrome, gamekeeper's thumb, tendonitis, thumb strains or sprains and post-cast removal
- > Constructed of vacuum-formed cross linked polyester with stretchable contact closure



PRODUCT NUMBER

	UNIV	PEDI	S/M	M/L
Universal	73-5000	73-5011		
Right			73-5001	73-5002
Left			73-5003	73-5004

SIZING GUIDE

UNIV	Fits most females & males
PEDI	Fits most children & young adolescents
S/M	Fits most females & smaller framed males
M/L	Fits most males & larger framed females

APPLICATION INSTRUCTIONS

1. Open straps.
2. Slide thumb into opening.
3. Secure straps around forearm.
4. Adjust thumb and forearm straps to desired comfort level.

CARE & CLEANING

Corflex manufactures a quality line of products and to help your care for the product properly we provided the following instructions for care & cleaning:

Remove any metal stays if possible, close all contact closures and hand wash in warm water with mild soap. Product should be air-dried, do not use high temperatures to wash or dry.

HOW TO PLACE AN ORDER

The Corflex Customer Service Team is ready to assist you in placing Your order. You may place your orders via a variety of methods:



TELEPHONE

1-800-426-7353, Ext1



FAX

1-603-623-4111



E-MAIL

Service@corflex.com



MAIL

CORFLEX, Inc.
Attn: Customer Service
669 East Industrial Park Dr
Manchester, NH 03109 USA

669 East Industrial Park Drive - Manchester, NH 03109
Phone (603) 623-3344 - Fax (603) 623-4111 - www.corflex.com

© Corflex | The advice and information contained in this document is not intended as a substitute for medical counseling. If you have any questions about treatment for a specific condition, please check with a health care professional.