



**CORFLEX**  
.....

**CORFLEX IS A GLOBALLY BASED  
MANUFACTURER AND DISTRIBUTOR OF QUALITY  
ORTHOPEDIC AND REHABILITATION PRODUCTS**



## ELBOW SLEEVES

### TARGET ELBOW SLEEVE

- > Provides support, warmth, protection and compression for the elbow
- > Indicated for treatment of elbow injuries and athletic protection
- > Manufactured from premium grade latex-free neoprene



### PRODUCT NUMBER

9"	S	M	L	XL	2XL
1/8"	88-3053	88-3054	88-3055	88-3056	88-3057
3/16"	88-3153	88-3154	88-3155	88-3156	88-3157

### SIZING GUIDE

S	M	L	XL	2XL
8"-9"	9"-10"	10"-11"	11"-12"	12"-14"

Measurement taken at widest part of forearm

### APPLICATION INSTRUCTIONS

1. With black outside exposed insert arm into the largest end of support and slide over elbow.
2. Make sure seam is on inside of arm and over bend of elbow.
3. Adjust placement as required.

### CARE & CLEANING

Corflex manufactures a quality line of products and to help your care for the product properly we provided the following instructions for care & cleaning:

Hand wash in warm water with mild soap. Product should be air-dried, do not use high temperatures to wash or dry.

### HOW TO PLACE AN ORDER

The Corflex Customer Service Team is ready to assist you in placing Your order. You may place your orders via a variety of methods:



**TELEPHONE**

1-800-426-7353, Ext1



**FAX**

1-603-623-4111



**E-MAIL**

Service@corflex.com



**MAIL**

CORFLEX, Inc.  
Attn: Customer Service  
669 East Industrial Park Dr  
Manchester, NH 03109 USA

669 East Industrial Park Drive - Manchester, NH 03109  
Phone (603) 623-3344 - Fax (603) 623-4111 - www.corflex.com

© Corflex | The advice and information contained in this document is not intended as a substitute for medical counseling. If you have any questions about treatment for a specific condition, please check with a health care professional.