



**CORFLEX**

**CORFLEX IS A GLOBALLY BASED  
MANUFACTURER AND DISTRIBUTOR OF QUALITY  
ORTHOPEDIC AND REHABILITATION PRODUCTS**



**PATELLA STABILIZER SUPPORTS / NEOPRENE**

**KNEE-O-TRAKKER**

- > Features textured C-shaped buttress for maximum patella control and versatile strapping system for dynamic treatment with adjustable tension control
- > Indicated for patella subluxation, patello-femoral syndrome, and as part of a rehab program for patella tendon injuries
- > Manufactured from premium grade latex-free neoprene with compression molded textured neoprene buttress



**PRODUCT NUMBER**

6"	XS	S	M	L	XL
3/16"	88-7232	88-7233	88-7234	88-7235	88-7236

6"	2XL	3XL	4XL
3/16"	88-7237	88-7238	88-7239

**SIZING GUIDE**

XS	S	M	L	XL
12"-13"	13"-14"	14"-15"	15"-16"	16"-18"

2XL	3XL	4XL
18"-20"	20"-22"	22"-24"

Measurement taken at mid-patella

**APPLICATION INSTRUCTIONS**

1. Loosen straps, making sure they remain threaded through d-rings.
2. After determining side of kneecap (patella) that requires support, position buttress on affected side and slide sleeve over knee.
3. Make sure seam is placed posteriorly in back of leg.
4. Your kneecap (patella) should be centered in the circular opening with buttress placed on rim of the kneecap (patella).
5. Grasp both straps by hook tip and pull back onto loop portion of strap.
6. Adjust straps to desired tension or comfort level.

**CARE & CLEANING**

Corflex manufactures a quality line of products and to help your care for the product properly we provided the following instructions for care & cleaning:

Close all contact closures and hand wash in warm water with mild soap. Product should be air-dried, do not use high temperatures to wash or dry.

**HOW TO PLACE AN ORDER**

The Corflex Customer Service Team is ready to assist you in placing Your order. You may place your orders via a variety of methods:



**TELEPHONE**

1-800-426-7353, Ext1



**FAX**

1-603-623-4111



**E-MAIL**

Service@corflex.com



**MAIL**

CORFLEX, Inc.  
Attn: Customer Service  
669 East Industrial Park Dr  
Manchester, NH 03109 USA

669 East Industrial Park Drive - Manchester, NH 03109  
Phone (603) 623-3344 - Fax (603) 623-4111 - www.corflex.com

© Corflex | The advice and information contained in this document is not intended as a substitute for medical counseling. If you have any questions about treatment for a specific condition, please check with a health care professional.